**Gel Impregnated**

**Ultramax<sup>TM</sup>**

*Knitted Double Velour **Vascular Graft***

**Instructions For Use**

**ATRIUM**

[www.atriummed.com](http://www.atriummed.com)

## **Description**

Ultramax™ is a knitted double velour vascular prosthesis made of Polyester yarn featuring an impregnation of meshed, absorbable gelatin of bovine origin. The gelatin is replaced by the ingrowth of fibrous tissue into the graft within 4 to 6 weeks.

## **Indications For Use**

Ultramax is indicated for use in reconstruction or replacement of damaged and diseased vessels of the abdomen in cases of aneurysmal or occlusive disease.

## **Contraindications**

None known.

## **Warnings And Precautions**

1. Caution: Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
2. Do not use after expiration date.
3. Do not resterilize Ultramax vascular prostheses.
4. Do not preclot. Ultramax vascular prostheses are impermeable to blood because of their impregnation.

## **Adverse Events**

Implantation of all vascular grafts have potential adverse patient effects or complications. These involve surgical risk from the operation, graft occlusion, coagulopathy, infection, thromboembolic episodes, aneurysm, hyperplasia, leaking or weeping of blood or serous fluid. In addition, coated grafts can cause immunological reactions and delayed healing. Any, and all of these complications can lead to amputation.

## **Clinical Studies**

Clinical data was collected from 95 patients with ectatic and obstructive disease of abdominal aorta in Japan. Postoperative complications associated with the diseases were found in 27% of the cases, (including: 3 cases each of pneumonia, hepatic dysfunction, renal dysfunction, bleeding, and infection; 1 case each of respiratory insufficiency, renal insufficiency, duodenal ulcer, bleeding, adhesive ileus, abdominal infection, abdominal hernia, thromboembolism of lower limb, deep vein embolism, swelling of right lower limb, thrombocytopenia, and stagnation of lymph fluid around graft.) None of these complications were deemed by the investigators to be graft related. Patency was satisfactory in all cases. Final observations at 15 months average (6 months to 25 months) found complications due to transplanted vascular prosthesis not present except one case, in which the graft on left leg was obstructed. In all the remaining cases, patency was confirmed. Therefore, cumulative patency rate was 98.8%.

## **Mode Of Application**

- The suitable shape, diameter and length of the prosthesis are selected according to the indication and anatomical conditions.
- The prosthesis must be sutured with round-bodied needles.
- If clamping of the prosthesis is necessary, the jaws of the clamp should be protected to avoid damaging the prosthesis.

## **Preclotting**

Ultramax vascular prosthesis are impermeable to blood because of their impregnation. They must not be preclotted.

## Sterilization And Storage

Ultramax vascular prosthesis are delivered gamma-sterilized. Ultramax grafts should be stored in their original packing at room temperature and a relative humidity of 50-60%.

## Resterilization

Ultramax vascular grafts must not be resterilized.

## Packaging

Ultramax vascular grafts come in a double-peel pack.

## SYMBOLS USED ON PRODUCT LABELS

**REF** CODE NUMBER    **LOT** LOT    **DIM** DIMENSIONS

**STERILE** **R** STERILE. STERILIZED BY GAMMA RADIATION.

 SEE PACKAGE INSERT     SINGLE USE ONLY

 EXPIRATION DATE    **Rx Only** PRESCRIPTION ONLY

 GEL IMPREGNATED

©2009

Ultramax and  
Atrium are trademark of  
Atrium Medical Corporation.

### ATRIUM MEDICAL CORPORATION

5 Wentworth Drive, Hudson, New Hampshire 03051 U.S.A.

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**ATRIUM EUROPE B.V.**
**ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.**