GB **Vascular Graft Tunneler**
Instructions for use

F

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ATRIUM

Product Description

The Atrium Tunneler is designed to create a subcutaneous tunnel for the placement of a vascular graft in both peripheral and dialysis applications. The Atrium Tunneler comes in three different rod configurations with a unique extender rod for extra long tunnels. Multiple diameter tunneler tips are available to create different sized tunnels. The Atrium Tunneler and its components are made of surgical-grade stainless steel and are designed for multiple reuse when properly maintained. If the threads of any component become damaged or do not fit properly, contact Atrium or its distributor to order replacement parts.

Indications For Use

The Atrium Tunneler is indicated for use as a device to create a subcutaneous tunnel for placement of a vascular prosthesis or autogenous graft for peripheral and extra-anatomic bypass procedures, and for vascular access.

Precautions

1. Avoid damaging the threaded portions of the rods and tips through improper handling.
2. The Atrium Tunneler requires disassembly, thorough cleaning, inspection, and sterilization prior to each use.

Tip Selection

To create the most appropriate size tunnel for the implant, use the following guidelines for tip selection. For regular, non-externally supported vascular grafts, natural vessels and similar products, use the same size tip as the internal diameter of the graft. For helix supported and ringed vascular grafts or similar products, use the next larger size tip, i.e., for a 6 mm Helix supported graft use a 7 mm tip. For tapered vascular grafts, consider the larger diameter end of the tapered graft for tip size selection.

Tunnel Rod Selection

The shaft of the standard Atrium Tunneler Rod is 6 mm in diameter. Selection of a particular tunneler rod model should be based on anatomical location, length and curve required. For a vascular tunnel longer than the standard tunneler length, use the Atrium extender rod by threading it directly into the tunneler shaft prior to tip attachment.

Operative Technique

For best results, always match the appropriate size tip to the graft diameter being inserted. Only hand tighten the threaded tip components. Once the tunneler has been passed into and through the tissue, and out through the exit incision, attach the graft over the tip and tie it securely with a suture. After graft attachment, pull the tunneler back through the tissue in a continuous and steady motion until the graft material exits the first incision. Avoid sudden stops and quick pulls as this may cause the graft to slip or tear off the tunneler tip. Cut the graft away from the tunneler and discard the portion of the material that was attached to the tunneler.

Cleaning And Sterilizing

The sterility of the product is the responsibility of the hospital. The following procedures are only recommendations and the sterility must be validated by the hospital.

For hospitals equipped with automatic washer/sterilizers:

Dissassemble all components of the tunneler and follow all instructions from the manufacturer of any automatic washer for correct handling and minimum exposure requirements for such surgical instruments.

For hospitals using steam sterilization:

Dissassemble all components of the tunneler and place all items in the designated location in the sterilization tray. Follow all instructions from the manufacturer of the steam sterilizer for correct handling and minimum exposure requirements for such surgical instruments.

For hospitals not equipped with a combination instrument and washer-sterilizer, the following cleaning instructions are recommended:

The Tunneler must be disassembled and rinsed in warm water. If the final wash is to be postponed, immerse in warm water solution containing a suitable detergent or blood solvent with a neutral pH (pH=7). DO NOT USE saline, bleach or harsh chemicals which may cause damage to the device. Use only a neutral pH (pH=7), low sudsing detergent for the final wash. DO NOT USE ordinary soap, which may leave an alkali film on the device. Ensure all adhering tissue and dried blood has been removed from all surfaces of the device including threads and handles. Inspect all components of the tunneler to ensure good working order and that they are clean prior to sterilization.

Replacement Parts / Ordering Information

Order No.	Description
26012	AV-1 Tunneler , 26 cm Length with 50° curve (Stainless Steel)
26013	AV-2 Tunneler , 27 cm Length with 180° curve (Stainless Steel)
26014	PV-1 Tunneler , 33 cm Length with 12° curve (Stainless Steel)
26011	XR-1 Tunneler Extender Rod , 25 cm Length (Stainless Steel)
26004	Sterilization Case (Polyphenylsulfone)
	Tunneler Tips (Stainless Steel)
26006	Set of 4 Size Tips (2 Per Size - 8 Total) 6 mm, 7 mm, 8 mm, 10 mm
26007	6 mm Tips (Set of 2)
26008	7 mm Tips (Set of 2)
26009	8 mm Tips (Set of 2)
26010	10 mm Tips (Set of 2)
26001	Complete Tunneler Set with Sterilization Case (3 Tunnelers, 1 Extender Rod, Set of 8 Tips)
26002	Complete Tunneler Set (3 Tunnelers, 1 Extender Rod, Set of 8 Tips)

SYMBOLS USED ON PRODUCT LABELS

REF CODE NUMBER **LOT** LOT NUMBER

 SEE PACKAGE INSERT **DIM** DIMENSIONS

Atrium is a trademark of Atrium Medical Corporation.
©2008 All Rights Reserved.

ATRIUM MEDICAL CORPORATION
5 Wentworth Drive, Hudson, New Hampshire 03051 U.S.A.
 603-880-1433  603-880-6718

ATRIUM EUROPE B.V.
Rendementsweg 20 B, 3641 SL Mijdrecht, The Netherlands
 +31-297-230-420  +31-297-282-653

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ATRIUM AUSTRALIA-PACIFIC RIM PTY.LTD.