# Pneumostat™ Chest Drain Valve

## Instructions For Use

USA

### Chest Drain Valve

GB

Valve de drainage thoracique

F

Thoraxdrainageventil

D

Válvula de drenaje torácico

E

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I

Válvula de drenagem torácica

P

Thoraxdrainklep

NL

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DK

Thoraxdräneringsventil

S

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FIN

Toraksdrenasjeventil

N

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
©2010  
Atrium and Pneumostat are trademarks  
of Atrium Medical Corporation  
Rev. 2010/04

001741

CE 0086

## Description

The Atrium Pneumostat™ is a disposable, single patient use chest drain valve with 30ml collection volume. A needleless fluid sample port provides fluid removal by syringe.

## Indications For Use

- To evacuate air from the chest cavity.

## Warnings

1. Do not use this device for fluid collection.
2. Do not obstruct air leak well.
3. Do not clamp patient tube during patient use.
4. Do not use this device for pneumonectomy.
5. Do not use or puncture needleless access port with needle.
6. Do not leave syringe attached to needleless fluid sample port.
7. Do not connect any open luer-lock connector to the needleless luer port located on the bottom of the chest drain valve.
8. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.

## Precautions

1. Replace Pneumostat chest drain valve when collected volume exceeds maximum capacity.
2. Replace Pneumostat chest drain valve if damaged.
3. Replace Pneumostat if evidence of occlusion has occurred.
4. Do not Y-connect two catheters to one Pneumostat.
5. Ensure Pneumostat is securely fastened to catheter.

## Set Up

Cut off beveled end of patient catheter and insert Pneumostat stepped connector end firmly into catheter. The Pneumostat stepped connector **E** fits thoracic catheter sizes 24Fr – 40Fr. The small stepped connector **F** and luer-lock connector **G** with tube are for smaller sized catheters.

## Air Leak Detection

Add 1ml water to air leak well to confirm air leak. Bubbling in water will confirm a patient air leak. Empty air leak well after use.

## Collection Graduations

Collection chamber graduations are calibrated in 5ml increments up to a maximum capacity of 30ml. Fluid graduations are accurate within 1ml with unit in a vertical position.

## Sampling Drainage

Must be in accordance with approved hospital infection control standards. Alcohol swab needleless fluid sample port prior to and after syringe attachment.

## System Disconnection

Clamp off all indwelling thoracic catheters prior to disconnecting chest drain valve from patient.

## System Disposal

Disposal of chest drain valve and contents should be in accordance with all applicable regulations.

## SYMBOLS USED ON PRODUCT PACKAGING

**REF** CODE NUMBER

**LOT** LOT NUMBER

**STERILE** **EO** STERILE. STERILIZED BY ETHYLENE OXIDE.

 SEE PACKAGE INSERT

 SINGLE USE ONLY

**Rx Only** PRESCRIPTION ONLY

 EXPIRATION DATE

Atrium and Pneumostat are registered trademarks of Atrium Medical Corporation.

©2010 All Rights Reserved.

CE 0086


**ATRIUM MEDICAL CORPORATION**  
5 Wentworth Drive, Hudson, New Hampshire 03051 U.S.A.  
 603-880-1433  603-880-6718

**ATRIUM EUROPE B.V.**  
Rendementsweg 20 B, 3641 SL Mijdrecht, The Netherlands  
 +31-297-230-420  +31-297-282-653

# Pneumostat™

## Chest Drain Valve

USA

### GB Features:

- A** One-Way Seal Valve
- B** Air Leak Well
- C** 30ml Collection Chamber
- D** Needleless Fluid Sample Port
- E** Stepped Connector
- F** Small stepped connector with tube
- G** Luer-lock connector with tube

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