
Gel Impregnated

*Ultramax*TM

*Knitted Double Velour **Vascular Patch***

Instructions For Use

ATRIUM

www.atriummed.com

Description

Ultramax™ Knitted Double Velour Vascular Patch is a warp-knitted double velour vascular patch made of polyester yarn featuring an impregnation of absorbable gelatin of bovine origin. The gelatin is replaced by the ingrowth of fibrous tissue in 4-6 weeks. The gelatin impregnation reduces the initial porosity of the patch from approximately 1900 ml/min/cm² to 0 ml/min/cm² at the time of implantation.

Indications

Ultramax Knitted Double Velour Vascular Patch is indicated for use in interventions in the regions of the carotid, profunda, femoral, and iliac arteries.

Warnings and Precautions:

1. Caution: Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
2. Do not use after expiration date.
3. Do not resterilize.
4. Do not preclot.

Adverse Events

Implantation of all vascular patches have potential adverse patient effects or complications. These involve surgical risk from the operation, vessel occlusion, coagulopathy, infection, thromboembolic episodes, aneurysm, hyperplasia, leaking or weeping of blood or serous fluid. In addition, coated patches can cause immunological reactions and delayed healing.

Mode of Action

The Ultramax Knitted Double Velour Vascular Patch is a vascular patch made of a porous polyester material which is impregnated with a biologically degradable coating (gelatin) to make the material blood tight.

This coating is enzymatically degraded and absorbed within 4-6 weeks after implantation. Endogenous tissue replaces the gelatin coating as it is absorbed and a cell layer, referred to as a neointima, forms on the inside of the patch. At the same time, connective tissue grows into the implant, anchoring it to the body of the recipient.

Contra-Indications

- Implantation in potentially infected areas
- Heart wall replacement
- Implantation in the region of the coronary arteries

Precautionary application

See mode of application.

Interactions with other agents

None known.

Mode of application

- The user must be familiar with the surgical technique, specific applications and the properties of the patch.
- Select the suitable Ultramax Knitted Double Velour Vascular Patch and cut the material to size in accordance with the indication and anatomical condition.
- Ultramax Knitted Double Velour Vascular Patch is not permeable to blood and should not be pre-clotted.
- Cut the patch to the required size using a sharp surgical instrument.

- Use only taper point suture needles as taper cut and other cutting needles may damage the patch material and/or host vessel.
- Clamping of the material should be avoided whenever possible and limited to clamps shod with soft material.

Side effects

None known.

Sterilization and Storage

Ultramax Knitted Double Velour Vascular Patches are delivered gamma-sterilized. They should be stored in their original packaging at room temperature and a relative humidity of 50-60%.

Resterilization

Ultramax Knitted Double Velour Vascular Patch must not be resterilized.

Packaging

Ultramax Knitted Double Velour Vascular Patch comes in a double-peel pack and should only be used if the packaging is undamaged.

SYMBOLS USED ON PRODUCT LABELS

REF CODE NUMBER **LOT** LOT NUMBER **DIM** DIMENSIONS

STERILE **R** STERILE. STERILIZED BY GAMMA RADIATION.

SEE PACKAGE INSERT

SINGLE USE ONLY

EXPIRATION DATE

Rx ONLY PRESCRIPTION ONLY

GEL IMPREGNATED

©2008

Ultramax and
Atrium are trademarks of
Atrium Medical Corporation.

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Hudson, New Hampshire 03051 U.S.A.
ATRIUM EUROPE B.V.
ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.