

> Nothing like it exists  
in the field of Hernia Repair


Clinical studies have shown  
that when it comes to  
Polypropylene Mesh

## “Less is More”


- ➔ 33% lighter with a smaller contact surface for minimal inflammation and an appropriate amount of connective tissue formation
- ➔ Superior plug stability helps prevent migration
- ➔ Innovative soft pile mesh enhances patient comfort

### Proven Design, ProLoop

The Most Biocompatible  
3 Dimensional Polypropylene  
Solution Available



*Increases tissue  
integration and  
promotes faster  
healing*




*Maximizes radial contact  
with surrounding tissue*

Lightweight or material reduced  
polypropylene presents significantly  
improved biocompatibility.



*Easy,  
self-centering  
insertion*



A Preferred Lighter Weight  
Fully Dimensional  
Hernia Plug



Klosterhalfen IL, Klinge U, Hermanns B, Schumpelick V. (2000) Pathology of traditional surgical nets for hernia repair after long-term implantation in humans. Der Chirurg 71:43-51



with Onlay >

Qty.	Part No.	Plug Height	Plug Width	Onlay Size	Box Qty.
—	30900	1.30" (3.3 cm)	1.10" (2.8 cm)	1.8 x 3.5" (4.6 x 8.9 cm)	2
—	30901	1.30" (3.3 cm)	1.55" (3.9 cm)	1.8 x 3.5" (4.6 x 8.9 cm)	2
—	30902	1.60" (4.1 cm)	1.90" (4.8 cm)	1.8 x 3.5" (4.6 x 8.9 cm)	2
—	30903	1.60" (4.1 cm)	2.0" (5.0 cm)	1.8 x 3.5" (4.6 x 8.9 cm)	2




with Keyhole Slit Onlay >

Qty.	Part No.	Plug Height	Plug Width	Onlay Size	Box Qty.
—	30908	1.30" (3.3 cm)	1.10" (2.8 cm)	2.0 x 3.5" (5.0 x 8.9 cm)	2
—	30909	1.30" (3.3 cm)	1.55" (3.9 cm)	2.0 x 3.5" (5.0 x 8.9 cm)	2
—	30910	1.60" (4.1 cm)	1.90" (4.8 cm)	2.0 x 3.5" (5.0 x 8.9 cm)	2
—	30911	1.60" (4.1 cm)	2.0" (5.0 cm)	2.0 x 3.5" (5.0 x 8.9 cm)	2



No Onlay >

Qty.	Part No.	Plug Height	Plug Width	Box Qty.
—	30920	1.30" (3.3 cm)	1.10" (2.8 cm)	2
—	30921	1.30" (3.3 cm)	1.55" (3.9 cm)	2
—	30922	1.60" (4.1 cm)	1.90" (4.8 cm)	2
—	30923	1.60" (4.1 cm)	2.0" (5.0 cm)	2



Signature \_\_\_\_\_ Date \_\_\_\_\_  
 For Evaluation Date \_\_\_\_\_  
 For Conversion Date \_\_\_\_\_



**ATRIUM MEDICAL CORPORATION**  
 5 Wentworth Drive  
 Hudson, New Hampshire 03051 U.S.A.  
 ☎ 603-880-1433 📠 603-880-6718

**ATRIUM EUROPE B.V.**  
 Rendementsweg 20B  
 3641 SL Mijdrecht, The Netherlands  
 ☎ +31-297-230-420 📠 +31-297-282-653

**ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.**  
 Level 6, 579 Harris Street  
 Ultimo NSW 2007 Australia  
 ☎ +61-2-8272-3100 📠 +61-2-8272-3199

© Atrium Medical Corporation 2010. All rights reserved. Printed in U.S.A. 10/10 Part #003356D. Atrium and ProLoop are trademarks of Atrium Medical Corporation. Marlex®, PerFix® is a registered trademark of DAVOL® Inc. Prolene® is a registered trademark of Ethicon® Inc.



Polypropylene Monofilament Mesh