


Vascular V12 covered stent



O R D E R I N G I N F O R M A T I O N

Standard Diameter Advanta™ V12 Covered Stent

Diameter	Length	Catheter Length	Introducer	Code #	<input type="checkbox"/>
5 mm	16 mm	80 cm	6 Fr	85340	<input type="checkbox"/>
5 mm	22 mm	80 cm	6 Fr	85341	<input type="checkbox"/>
5 mm	38 mm	80 cm	7 Fr	85320	<input type="checkbox"/>
5 mm	59 mm	80 cm	7 Fr	85321	<input type="checkbox"/>
6 mm	16 mm	80 cm	6 Fr	85342	<input type="checkbox"/>
6 mm	22 mm	80 cm	6 Fr	85343	<input type="checkbox"/>
6 mm	38 mm	80 cm	7 Fr	85322	<input type="checkbox"/>
6 mm	59 mm	80 cm	7 Fr	85323	<input type="checkbox"/>
7 mm	16 mm	80 cm	7 Fr	85344	<input type="checkbox"/>
7 mm	22 mm	80 cm	7 Fr	85345	<input type="checkbox"/>
7 mm	38 mm	80 cm	7 Fr	85324	<input type="checkbox"/>
7 mm	59 mm	80 cm	7 Fr	85325	<input type="checkbox"/>
8 mm	38 mm	80 cm	7 Fr	85326	<input type="checkbox"/>
8 mm	59 mm	80 cm	7 Fr	85327	<input type="checkbox"/>
9 mm	38 mm	80 cm	7 Fr	85328	<input type="checkbox"/>
9 mm	59 mm	80 cm	7 Fr	85329	<input type="checkbox"/>
10 mm	38 mm	80 cm	7 Fr	85360	<input type="checkbox"/>
10 mm	59 mm	80 cm	7 Fr	85361	<input type="checkbox"/>

Diameter	Length	Catheter Length	Introducer	Code #	<input type="checkbox"/>
5 mm	16 mm	120 cm	6 Fr	85350	<input type="checkbox"/>
5 mm	22 mm	120 cm	6 Fr	85351	<input type="checkbox"/>
5 mm	38 mm	120 cm	7 Fr	85330	<input type="checkbox"/>
5 mm	59 mm	120 cm	7 Fr	85331	<input type="checkbox"/>
6 mm	16 mm	120 cm	6 Fr	85352	<input type="checkbox"/>
6 mm	22 mm	120 cm	6 Fr	85353	<input type="checkbox"/>
6 mm	38 mm	120 cm	7 Fr	85332	<input type="checkbox"/>
6 mm	59 mm	120 cm	7 Fr	85333	<input type="checkbox"/>
7 mm	16 mm	120 cm	7 Fr	85354	<input type="checkbox"/>
7 mm	22 mm	120 cm	7 Fr	85355	<input type="checkbox"/>
7 mm	38 mm	120 cm	7 Fr	85334	<input type="checkbox"/>
7 mm	59 mm	120 cm	7 Fr	85335	<input type="checkbox"/>
8 mm	38 mm	120 cm	7 Fr	85336	<input type="checkbox"/>
8 mm	59 mm	120 cm	7 Fr	85337	<input type="checkbox"/>
9 mm	38 mm	120 cm	7 Fr	85338	<input type="checkbox"/>
9 mm	59 mm	120 cm	7 Fr	85339	<input type="checkbox"/>
10 mm	38 mm	120 cm	7 Fr	85364	<input type="checkbox"/>
10 mm	59 mm	120 cm	7 Fr	85365	<input type="checkbox"/>

NEW > Large Diameter Advanta™ V12 Covered Stent

Diameter	Length	Catheter Length	Introducer	Code #	<input type="checkbox"/>
12 mm	29 mm	80 cm	9 Fr	85370	<input type="checkbox"/>
12 mm	41 mm	80 cm	9 Fr	85371	<input type="checkbox"/>
12 mm	61 mm	80 cm	9 Fr	85372	<input type="checkbox"/>
*14 mm	29 mm	80 cm	11 Fr	85373	<input type="checkbox"/>
*14 mm	41 mm	80 cm	11 Fr	85374	<input type="checkbox"/>
*14 mm	61 mm	80 cm	11 Fr	85375	<input type="checkbox"/>
*16 mm	29 mm	80 cm	11 Fr	85376	<input type="checkbox"/>
*16 mm	41 mm	80 cm	11 Fr	85377	<input type="checkbox"/>
*16 mm	61 mm	80 cm	11 Fr	85378	<input type="checkbox"/>

Diameter	Length	Catheter Length	Introducer	Code #	<input type="checkbox"/>
12 mm	29 mm	120 cm	9 Fr	85379	<input type="checkbox"/>
12 mm	41 mm	120 cm	9 Fr	85380	<input type="checkbox"/>
12 mm	61 mm	120 cm	9 Fr	85381	<input type="checkbox"/>
*14 mm	29 mm	120 cm	11 Fr	85382	<input type="checkbox"/>
*14 mm	41 mm	120 cm	11 Fr	85383	<input type="checkbox"/>
*14 mm	61 mm	120 cm	11 Fr	85384	<input type="checkbox"/>
*16 mm	29 mm	120 cm	11 Fr	85385	<input type="checkbox"/>
*16 mm	41 mm	120 cm	11 Fr	85386	<input type="checkbox"/>
*16 mm	61 mm	120 cm	11 Fr	85387	<input type="checkbox"/>

*Pending CE Approval

Notes: _____

Signature: _____ Date: _____



ATRIUM MEDICAL CORPORATION
5 Wentworth Drive
Hudson, New Hampshire 03051 U.S.A.
☎ 603-880-1433 📠 603-880-6718

ATRIUM EUROPE B.V.
Rendementsweg 20 B
3641 SL Mijdrecht, The Netherlands
☎ +31-297-230-420 📠 +31-297-282-653

ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.
L 1 Bridgepoint, 3 Brady Street
Mosman NSW 2088 Australia
☎ +61-2-9960-0169 📠 +61-2-8969-2735