




## The No Weep Guarantee™

### > An Industry First


The all new **FLIXENE** “next generation” composite vascular graft. Utilizing Atrium’s state-of-the-art Trilaminar Composite Construction (TCC), **FLIXENE** has been engineered to surpass the toughest clinical demands required in today’s more challenging vascular procedures. Whether for extreme dialysis access conditions or lower limb reconstructions, **FLIXENE** provides unmatched strength, durability, weep protection and security.

**FLIXENE**’s exclusive hydrostatic protection membrane minimizes hydrostatic pressure, thereby eliminating blood strikethrough. This hydrostatic protection membrane provides **FLIXENE** with the highest water entry pressure of any vascular graft commercially available. The result is a product that is the least likely of any vascular graft to experience weeping through the graft wall after implantation.

### Trilaminar Composite Construction



### Impact of Cannulation on Radial Tensile Strength of Vascular Grafts\*



**FLIXENE**’s unique 3-layer construction provides the best “mechanical solution” for your most challenging vascular repairs. Atrium’s proprietary structure is designed to withstand the rigors of multiple needle cannulations, kink and torque compression resistance, aneurysm formation and/or diameter creep following implantation. At installation, **FLIXENE** is twice as strong as the nearest competitor. After 20 needlesticks per cm with a 14 gauge needle, testing demonstrates **FLIXENE** is 3 times stronger than most PTFE products. **FLIXENE** is designed to take the wear and tear that today’s vascular grafts require from dialysis access to lower limb reconstruction.

**FLIXENE** is one of the the most durable vascular implants available. By implanting **FLIXENE**, you will have the confidence that you are providing the best surgical intervention solution for your patient.

### Guarantee Requirements

Atrium is so confident in your experience with the **FLIXENE** vascular graft that it is providing a NO WEEP SATISFACTION GUARANTEE. It is guaranteeing your satisfaction that the product will not weep when used as directed. If the product does weep, Atrium will provide the hospital with a no charge replacement.

### Conditions:

The only conditions are the following:

- 1) Fluid must not be forced through the graft wall using a syringe or other method.
- 2) The physician must release the distal clamp first to “back bleed” the graft then release the proximal clamp. This prevents the “hammer effect” which is a common cause of weeping in vascular grafts.
- 3) The graft must not be pre-wet by exposure to organic solvents such as alcohol.

### Fulfillment of the Guarantee

The implanting physician must complete and sign the **FLIXENE** Guarantee form. Once the form is completed and provided to the Atrium territory representative, Atrium will provide one free replacement graft, not cash or credit to the hospital, for the physician’s next case after receipt of the corresponding implant tracking card. Replacement product will only be applicable for those graft items that have been purchased by the hospital.



# Durability That Instills Confidence