










Chest Drain Accessories

Item #	Description
19910 	Red Cap <ul style="list-style-type: none"> • For capping patient line before disposal • For use if only one patient line is needed on a dual drain • Twenty-five (25) each per case – sterile
19913 	Adult/Pediatric Connector 1/4" x 3/8" <ul style="list-style-type: none"> • For use with an adult drain and smaller catheters – 16, 20 French • Twenty-five (25) each per case – sterile
19914 	Standard Adult Connector 3/8" x 3/8" <ul style="list-style-type: none"> • Standard 5-in-1 connector • For use on larger catheters – 24 French and larger • Twenty-five (25) each per case – sterile
19915 	Pediatric Connector 1/4" x 1/4" <ul style="list-style-type: none"> • Larger connector included in Pedi-drain • For use with an adult drain and smaller catheters – 16 French or larger • Twenty-five (25) each per case – sterile
19916 	Sims Connector 3/16" x 1/4" <ul style="list-style-type: none"> • Same connector included in Pedi-drain • For use with smaller catheter – 8, 10, 12 French • Twenty-five (25) each per case – sterile
19917 	Y-Connector 3/8" x 3/8" x 3/8" <ul style="list-style-type: none"> • Connects 2 chest tubes to one adult chest drain • Twenty-five (25) each per case – sterile
19918 	Step-Down Connector 3/16" x 3/8" <ul style="list-style-type: none"> • For use with adult chest drain • For use with Blake Drain to adult drain • Fits smaller catheters – 8, 10, 12 French • Twenty-five (25) each per case – sterile
19907 	Latex-Free Tube Pack <ul style="list-style-type: none"> • Safety stock if patient tube sterility is accidentally compromised • Provides cost savings to patient and hospital • Converts a latex drain for latex sensitive patients • 3/8" ID tubing • Ten (10) each per case - sterile
19908 	Latex Tube Pack <ul style="list-style-type: none"> • Safety stock if patient tube sterility is accidentally compromised • Provides cost savings to patient and hospital • 3/8" ID tubing • Ten (10) each per case – sterile

Chest Drain Accessories

19909

Y-Extension Sets



- Simple, pre-assembled extension connection for attaching 2 chest tubes to one chest drain
- Ten (10) each per case – sterile
- Latex-free

19921

Needleless Access Sampling Port



- Allows for fast, safe specimen collection
- Attaches to inline connector
- Latex-free
- Ten (10) each per case – sterile

19927

Luer-Lock Connector With Tube



- Allows connection of a standard chest drain or mobile drain to a luer-lock pneumothorax catheter
- Latex-free
- Twenty-five (25) each per case – sterile

19928

Suction Y-Adaptor



- For use when connecting two chest drains to one suction source
- Latex-free
- Ten (10) each per case – sterile

19929

PleuraGuide™ Chest Tube Introducer



- For use with catheters 24 French or larger
- Allows for easier, less traumatic chest tube insertion
- Latex-free
- Twenty (20) each per case – sterile

19930

Express Suction Bulb



- For use with the Express Mini 500
- Allows portable, manual vacuum application
- Latex-free
- Ten (10) each per case – sterile

19911 Standard

Drain-Caddy



- Eliminates knockover
- Quick, easy, secure attachment
- Long lasting, reusable
- Use multiple units on one IV pole
- One (1) each per case

Atrium Medical Corporation
5 Wentworth Drive
Hudson, NH 03051 USA
Phone 603-880-1433
Fax 603-880-6718

Atrium Europe B.V.
Rendementsweg 20 B
3641 SL Mijdrecht, The Netherlands
Phone +31 297 230 420
Fax +31 297 282 653

Atrium Australia-Pacific Rim Pty. Ltd.
Level 6, 579 Harris Street
Ultimo NSW 2007
Phone +61 8272 3100
Fax +61 8272 3199



www.atriummed.com