



Join the Revolution.

Express™

DRY SEAL CHEST DRAIN

To help chest tube patients get ambulatory as soon as possible, Atrium is proud to offer the **Express™ Mini 500 Dry Seal Chest Drain**. The Express Mini 500 is a waterless operating system with 500ml collection volume, a dry suction regulator preset at $-20\text{cmH}_2\text{O}$, and a dry one-way valve for seal protection.

The Express Mini 500 is the world's 1st fully-wearable chest drain to get your patients mobile faster and safer.

Atrium is leading the revolution in mobile chest drainage. Studies have shown there are advantages to getting patients with chest tubes ambulatory as soon as possible. Early ambulation can also play a role in reducing length of hospital stay, which equates to dollars saved by the hospital. Are you ready to get mobile?

Simple Air Leak Detection

MINI
500



Vacuum Indicator

500ml Collection Capacity

Empty Collection Chamber
Via Needleless Luer Port

Ordering Information

Part Number	Description	Case Quantity
16400	Express Mini 500	6
19930	Express Suction Bulb	10



ATRIUM MEDICAL CORPORATION
5 Wentworth Drive
Hudson, New Hampshire 03051 U.S.A.
☎ 603-880-1433 ☎ 603-880-6718

ATRIUM EUROPE B.V.
Rendementsweg 24 L
3641 SL Mijdrecht, The Netherlands
☎ +31-297-230-420 ☎ +31-297-282-653

ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.
L 1 Bridgepoint, 3 Brady Street
Mosman NSW 2088 Australia
☎ +61-2-9960-0169 ☎ +61-2-8969-2735