


Firm PVC Thoracic Catheters

Smooth, atraumatic tip
allows faster, easier insertion.



Smooth
**ATRAUMATIC
TIP**



Atrium Firm PVC Thoracic Catheters

- ▶ Open tip provides easy stab wound placement and fast, secure connection to chest drain tubing.
- ▶ Kink-resistant, thermosensitive PVC conforms easily to anatomical contour for superior patient comfort and improved drainage performance.
- ▶ 2 cm depth markings aid in accurate catheter placement.
- ▶ Gentle spiral hole configuration helps reduce tissue invagination for optimal drainage efficiency.
- ▶ Tapered oval drainage eyes maximize drainage performance and help prevent tissue damage and clot adhesion.
- ▶ Blue radiopaque stripe confirms placement upon X-ray.
- ▶ Complete line of products.



Advanced Extrusion Design Improves Patient Comfort, Surgical Placement, and Provides a Preferred Atraumatic Insertion Technique



- ▶ **Smooth, Atraumatic Tip Design**
 - Allows smooth, easy insertion with less trauma
 - Full depth markings for accurate placement
 - Convenient to grab, open tip
 - Provides fast, secure catheter connection

- ▶ **Tapered Oval Drainage Eyes**
 - Gentle, spiral hole configuration helps reduce tissue invagination
 - Helps reduce mechanical irritation, tissue damage, and clot adhesion
 - Provides optimal drainage efficiency

Ordering Information

Firm PVC Catheters

Description	Length	# Eyes	Code #	UPN #
12FR Straight	14"	4	15012	20650862150128
16FR Straight	14"	4	15016	20650862150166
20FR Straight	19"	4	15020	20650862150203
24FR Straight	19"	5	15024	20650862150241
28FR Straight	20"	6	15028	20650862150289
32FR Straight	20"	6	15032	20650862150326
36FR Straight	20"	6	15036	20650862150364
40FR Straight	20"	6	15040	20650862150401
12FR Right Angle	14"	4	15212	20650862152122
16FR Right Angle	14"	4	15216	20650862152160
20FR Right Angle	19"	4	15220	20650862152207
24FR Right Angle	19"	4	15224	20650862152245
28FR Right Angle	20"	5	15228	20650862152283
32FR Right Angle	20"	5	15232	20650862152320
36FR Right Angle	20"	5	15236	20650862152368
40FR Right Angle	20"	5	15240	20650862152405

Atrium Firm PVC Catheters are packaged sterile 1 per pouch, 10 pouches per case.

Sterile • Non-Pyrogenic

Federal (USA) law restricts this device to sale by or on the order of a physician.

Visit us on the web at: www.atriummed.com

Atrium Medical Corporation 5 Wentworth Drive, Hudson, New Hampshire U.S.A. 03051
Telephone 603-880-1433 1-800-528-7486 FAX 603-880-6718

European Representative:
Rendementsweg 24 L, 3641 Mijdrecht, The Netherlands

☎ 31 297 230 420 ☒ 31 297 230 422

