


NEW!

HydraGlide™ XL

SILICONE THORACIC CATHETERS

Product Features

- Made of soft, kink resistant silicone.
- Extended length catheter with 15-18 eyelets for maximum drainage efficiency.
- Smaller bore 20Fr and 24Fr sizes.
- Proprietary HydraGlide™ coating features a covalently bonded Heparin complex which reduces adherent clot formation in and around the catheter for superior drainage performance.
- Hydrophilic surface coating with its low coefficient of friction enhances patient comfort, particularly upon tube removal.
- Designed specifically for cardiac and thoracic drainage applications.
- Proven single lumen design drains more than twice the amount of fluid compared to similar size multi-channel wound drains.
- High tear strength silicone reduces risk of tearing or breaking.
- Croppable to reduce length, if required.
- Blue radiopaque stripe confirms placement upon x-ray.



**To place an order,
call 1-800-528-7486.**



www.atriummed.com

© Atrium Medical Corporation 2003.
Made and printed in U.S.A. 12/03 003307