


HydraGlide™ Coated Silicone Chest Tube Catheters



Flat
Mediastinal
Chest Tube
Sizes

NEW
SIZES

For enhanced drainage and patient comfort...



Atrium **HydraGlide™** Silicone Chest Tube and Mediastinal Catheters

- ▶ Permanent, covalently bonded heparin complex inhibits clot formation.
- ▶ Ultra low coefficient of friction helps reduce mechanical irritation to heart and lung tissues.
- ▶ Superior strength reduces risk of tearing or breaking.
- ▶ Clear, premium quality silicone provides enhanced drainage visibility.
- ▶ Soft, kink-resistant silicone conforms easily to anatomy for improved patient comfort.
- ▶ Blue radiopaque stripe confirms catheter placement upon X-ray.



Ordering Information

New!

HydraGlide™ XL Silicone Chest Tube Catheters

Description	Length	# Eyes	Order #
20FR Extended Length - Straight	26"	15	14720
24FR Extended Length - Straight	26"	18	14724

HydraGlide™ Silicone Chest Tube Catheters

Description	Length	# Eyes	Order #
16FR Straight	20"	4	14116
20FR Straight	20"	5	14120
24FR Straight	20"	6	14124
28FR Straight	20"	6	14128
32FR Straight	20"	6	14132
36FR Straight	20"	6	14136

HydraGlide™ Silicone Mediastinal Catheters

Description	Length	# Eyes	Order #
7mm Flat, Eyes on Top & Bottom	20"	6	14507
9mm Flat, Eyes on Top & Bottom	20"	9	14509
7mm Flat, Eyes on Top	20"	4	14707
9mm Flat, Eyes on Top	20"	6	14709
7mm Flat, Eyes on Side	20"	5	14907
9mm Flat, Eyes on Side	20"	5	14909

Atrium HydraGlide™ silicone chest tube and mediastinal catheters are packaged sterile 1 per pouch, 10 pouches per case.

Sterile • Non-Pyrogenic

Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING

This product contains a permanent, covalently bonded Heparin complex to inhibit clot formation. It is now known that selected patients may be sensitive to Heparin exposure following prior Heparin usage. Patients must be screened for Heparin induced thrombocytopenia (HIT) prior to use of this device. **Do not use this device if any patient is known to be sensitive to Heparin or is suspected of HIT at any time.** Please refer to Atrium's website for additional references and information regarding HIT.



Atrium Medical Corporation
5 Wentworth Drive
Hudson, NH 03051 USA
Phone 603-880-1433
Fax 603-880-6718

Atrium Medical Corporation
Remdementsweg 20 B
3641 SL Mijdrecht, The Netherlands
Phone +31 297 230 420
Fax +31 297 282 653

Atrium Australia-Pacific Rim Pty. Ltd.
L1 Bridgepoint
3 Brady Street
Mosman NSW 2088 Australia
Phone +61 2 9960 0169
Fax +61 2 8969 2735