

# Oasis™

## DRY SUCTION CHEST DRAINS

**The most  
Advanced  
water seal  
chest drain  
Available.**

Atrium has once again proven itself as a true innovator by fulfilling previously unmet clinical needs in a way that benefits both the clinician and patient. By including pre-packaged sterile water to fill the water seal along with unsurpassed knock-over protection, the Oasis dry suction chest drain offers the ultimate in fast setup and patient safety.

The next generation Oasis combines the safe, dependable vacuum control of its state-of-the-art dry suction regulator with a fully calibrated water seal. Other features include multi-position hangers, latex-free tubing and 24 hour technical support...because emergency questions do not work an 8 hour day.



**PRE-PACKAGED\*  
WATER**

**In-Line Connector  
with needless  
access sampling port**

**New Larger Easy-  
to-Grab Handle  
and Flexi-Hangers  
allow the drain to be  
hung bedside from a  
single point**


**Convenient  
Fingertip  
Suction  
Adjustment**

**Pre-Packaged  
Water  
offers the ultimate in  
fast, convenient setup**

**Innovative Knock-  
Over Nozzles  
allow fluid levels to be  
recovered in the event  
of an accidental knock-  
over**

**Large, Easy-to-  
Read Graphics  
provide fast,  
accurate drainage  
assessment**

**New Graduated Air  
Leak Monitor  
for assessing patient  
air leak trends**



\*Patented

  
**SPEED. RECOVERY.**

 **ATRIUM**

# Oasis™

## DRY SUCTION CHEST DRAINS

### Atrium® Oasis™ Dry Suction Chest Drains

#### Ordering Information

Model	Description	Case Quantity
<b>3600-100</b>	Single Collection, 1 Patient Tube ATS Bag Compatible	6
<b>3600-150</b>	Single Collection, 1 Patient Tube Pre-connected ATS Bag	3
<b>3620-100</b>	Dual Collection, 2 Patient Tubes ATS Bag Compatible,	6
<b>3612-400</b>	Infant/Pediatric, 1 Patient Tube 1/4" Tubing, Pediatric Connectors	6
<b>3650-100</b>	ATS Blood Recovery, 1 Patient Tube ATS Bag Compatible	6
<b>2550</b>	In-Line ATS Blood Bag For use with Atrium ATS Compatible Chest Drains	6
<b>2450</b>	Self-Filling ATS Blood Bag For use with Atrium ATS Blood Recovery Drains	10
<b>19911</b>	Drain-Caddy™ Chest Drain I.V. Pole Adapter	1

Federal (USA) law restricts this device to sale by or on the order of a physician. Please refer to our detailed product catalog for a complete listing of models and pricing. To place an order or to request product information call Atrium toll free at 1-800-528-7486.



#### Flexible Options for Post-Op ATS

Whether for routine ATS or for emergency chest trauma applications, Atrium's ATS Blood Recovery models offer unlimited flexibility for prescribing continuous pump reinfusion or post-event blood retrieval with a 2450 Self-Filling ATS blood bag. Models equipped with in-line connectors allow rapid in-line ATS blood bag attachment when required.



**ATRIUM MEDICAL CORPORATION**  
5 Wentworth Drive  
Hudson, New Hampshire 03051 U.S.A.  
☎ 603-880-1433 📠 603-880-6718

**ATRIUM EUROPE B.V.**  
Rendementsweg 20 B  
3641 SL Mijdrecht, The Netherlands  
☎ +31-297-230-420 📠 +31-297-282-653

**ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.**  
L 1 Bridgepoint, 3 Brady Street  
Mosman NSW 2088 Australia  
☎ +61-2-9960-0169 📠 +61-2-8969-2735