

Pneumostat™
Chest Drain Valve.
The Future of Mobile
Patient Recovery.



Studies have shown there are many advantages associated with getting patients with chest tubes ambulatory as soon as possible. These include: enhanced circulation, particularly to the lower extremities; earlier restoration of normal pulmonary function; improved patient mood with increasing independence; and decreased need for sedation and pain control drugs. Early ambulation can also play a role in reducing length of hospital stay, which equates to dollars saved by the hospital.

That's why Atrium developed the Pneumostat Chest Drain Valve. Pneumostat is a total package, unique from other chest drain valves. Whether you have the need to treat a patient with simple pneumothorax or you are transitioning a patient from a traditional chest drain to a smaller, more mobile device, Pneumostat gives you and your patients *the ultimate low cost solution.*

The
Ultimate
Freedom
Machine



Pneumostat™
Chest Drain Valve

Ordering Information

Description	Code	UPN Number	Qty.
Pneumostat™ Chest Drain Valve w/ 30ml Collection	16100	20650862161001	10/Bx

To place an order or request product information, call Atrium toll free at **1-800-528-7486.**



Visit us on the web at:

www.atriummed.com

ATRIUM MEDICAL CORPORATION

5 Wentworth Drive
Hudson, New Hampshire 03051 U.S.A.
☎ 603-880-1433 📠 603-880-6718

ATRIUM EUROPE

Rendementsweg 24, L III
3641 SL Mijdrecht, The Netherlands
☎ +31-297-230-420 📠 +31-297-230-422

ATRIUM AUSTRALIA-PACIFIC RIM PTY.LTD.

L 1 Bridgepoint
3 Brady Street
Mosman NSW 2088 Australia
☎ +61-2-9960-0169 📠 +61-2-8969-2735



Pneumostat™
Chest Drain Valve

Get Your Patients Mobile!



Freedom
Machines



**Fast, secure
garment
attachment**

**Secure
connection
to any size
thoracic catheter**

**Ready to
use, no
priming
needed**

**State-of-the-art
dry seal valve
provides enhanced
patient protection**

**Unique Air
Leak Well
provides
quick, simple
air leak
confirmation**

**Lightweight,
low-profile
design enhances
patient comfort
and mobility**

**Closed,
easy to
empty 30ml
collection
chamber**

**Luer port allows
safe, convenient
fluid removal**

**100%
Latex
Free**

Pneumostat™

Chest Drain Valve

- ▶ **Fast, Easy Attachment to Any Size Catheter**
- ▶ **Closed System is Easy to Empty**
- ▶ **Simple Air Leak Detection**
- ▶ **Needleless Luer Port for Convenient Fluid Management**

 **ATRIUM**