

Thermosensitive PVC Chest Tube Catheters



Taper tip design allows faster,
easier atraumatic insertion.



Atrium *Taper Tip* Chest Tube Catheters

- ▶ Advanced tapered tip provides easy stab wound placement and fast, secure connection to chest drain tubing.
- ▶ Soft, kink-resistant, thermosensitive PVC conforms easily to anatomical contour for superior patient comfort and improved drainage performance.
- ▶ 2 cm depth markings aid in accurate catheter placement.
- ▶ Gentle spiral hole configuration helps reduce tissue invagination for optimal drainage efficiency.
- ▶ Tapered oval drainage eyes maximize drainage performance and help prevent tissue damage and clot adhesion.
- ▶ Blue radiopaque stripe confirms placement upon X-ray.
- ▶ Ideal for all ATS considerations.

100%
Latex Free



Advanced Extrusion Design Improves Patient Comfort, Surgical Placement, and Provides a Preferred Atraumatic Insertion Technique



▶ Atraumatic Taper Tip Design

- Allows smooth, easy insertion with less trauma
- Provides clean, safe, contaminant-free lumen and catheter connection
- Full depth markings for accurate placement
- Convenient to grab, open tip

▶ ATP Trocar Catheters

- Tapered ATP guiding stylet
- Convenient catheter placement
- Available in adult, pediatric, and infant sizes

▶ Tapered Oval Drainage Eyes

- Gentle, spiral hole configuration helps reduce tissue invagination
- Helps reduce mechanical irritation, tissue damage, and clot adhesion
- Provides optimal drainage efficiency

Ordering Information

PVC Catheters

Description	Length	# Eyes	Code #
8FR Straight	8"	2	8008
10FR Straight	8"	2	8010
12FR Straight	16"	4	8012
16FR Straight	16"	4	8016
20FR Straight	22"	4	8020
24FR Straight	22"	5	8024
28FR Straight	23"	6	8028
32FR Straight	23"	6	8032
36FR Straight	23"	6	8036
40FR Straight	23"	6	8040
12FR Right Angle	16"	4	8112
16FR Right Angle	16"	4	8116
20FR Right Angle	22"	4	8120
24FR Right Angle	22"	4	8124
28FR Right Angle	23"	5	8128
32FR Right Angle	23"	5	8132
36FR Right Angle	23"	5	8136
40FR Right Angle	23"	5	8140
8FR Trocar	8"	2	8408
10FR Trocar	8"	2	8410
12FR Trocar	9"	2	8412
16FR Trocar	10"	2	8416
20FR Trocar	16"	2	8420
24FR Trocar	16"	2	8424
28FR Trocar	16"	2	8428
32FR Trocar	16"	2	8432

Atrium Thermosensitive PVC Catheters are packaged sterile 1 per pouch, 10 pouches per case.

Sterile • Non-Pyrogenic

Federal (USA) law restricts this device to sale by or on the order of a physician.



Atrium Medical Corporation
5 Wentworth Drive
Hudson, NH 03051 USA
Phone 603-880-1433
Fax 603-880-6718

Atrium Medical Corporation
Rendementsweg 24, L111
3641 SL Mijdrecht, The Netherlands
Phone +31 297 230 420
Fax +31 297 230 422

Atrium Australia-Pacific Rim Pty. Ltd.
L1 Bridgepoint
3 Brady Street
Mosman NSW 2088 Australia
Phone +61 2 9960 0169
Fax +61 2 8969 2735

Visit Our Website:

www.atriummed.com