

*Chest Drains, Thoracic Catheters,  
Chest Drain Valves and Accessories*



**ATRIUM**

**Product Catalog**

Atrium Medical Corporation  
Effective April 2008

**Atrium Medical Corporation**

5 Wentworth Drive  
Hudson, NH 03051 USA

☎ (603)880-1433

☎ (603)880-6718

**Atrium Europe B.V.**

Rendementsweg 20 B  
3641 SL Mijdrecht  
The Netherlands

☎ +31 297 230 420

☎ +31-297-282-653

**Atrium Australia-Pacific Rim Pty. Ltd.**

Sales Office  
L1 Bridgepoint  
3 Brady Street  
Mosman NSW 2088

☎ (02) 9960-0169

☎ (02) 8969-2735

*www.atriummed.com*

## Ocean Water Seal Chest Drains

### Single Collection

♦ 1 Patient Tube, Sterile

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
2002-000	A7000L, A7000LF A7000-08L, A7000-08LF	8888571315	8884713176	6	<b>Single Collection - ATS Compatible</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock
2002-000L	A7000LF	571315	713176	6	<b>Single Collection - ATS Compatible</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock ♦ Latex Patient Tube
2002-070	N/A	N/A	N/A	6	<b>Single Collection - ATS Compatible</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock ♦ Sterile Fluid Path Packaging
2002-100	A7000L, A7000LF A7000-08L, A7000-08LF	8888571315	8884713176	6	<b>Single Collection - ATS Compatible</b> ♦ In-Line Patient Tube Connector & Clamp ♦ No Suction Control Stopcock
2002-300	A8000L, A8000LF A8000-08L, A8000-08LF	8888571299	8884713308	6	<b>Single Collection - Standard System</b> ♦ No In-Line Patient Tube Connector ♦ Suction Control Stopcock
2002-400	A8000L, A8000LF A8000-08L, A8000-08LF	8888571299	8884713308	6	<b>Single Collection - Standard System</b> ♦ No In-Line Patient Tube Connector ♦ No Stopcock ♦ Pre-packaged Sterile Water ♦ No Suction Tubing
2002-040	N/A	N/A	N/A	6	<b>Single Collection - ATS Compatible</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock ♦ 1/4" x 1/4" Pedi-Connector
2002-050	A7050	8888571323	8888713192	3	<b>Single Collection - ATS Ready</b> ♦ Preattached In-Line ATS Blood Bag ♦ Suction Control Stopcock



## Ocean Water Seal Chest Drains

### Single Collection

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
2002-050L	A7050	571323	713192	3	<b>Single Collection - ATS Ready</b> ♦ Preattached In-Line ATS Blood Bag ♦ Suction Control Stopcock ♦ Latex Patient Tube
2002-057	N/A	N/A	N/A	3	<b>Single Collection - ATS Ready</b> ♦ Preattached In-Line ATS Blood Bag ♦ Suction Control Stopcock ♦ Sterile Fluid Path Packaging

### Dual Collection

♦ 2 Patient Tubes, Major/Minor Chambers, Sterile

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
2020-000	A7002L, A7002LF A7002-08L, A7002-08LF	N/A	N/A	6	<b>Dual Collection - ATS Compatible</b> ♦ 2 In-Line Patient Tube Connectors & Clamps ♦ Suction Control Stopcock
2020-100	A7002L, A7002LF A7002-08L, A7002-08LF	N/A	N/A	6	<b>Dual Collection - ATS Compatible</b> ♦ 2 In-Line Patient Tube Connectors & Clamps ♦ No Suction Control Stopcock
2020-300	A8002L, A8002LF A8002-08L, A8002-08LF	8888571406	N/A	6	<b>Dual Collection - Standard System</b> ♦ No In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock



## Ocean Water Seal Chest Drains

### Infant/Pediatric

♦ 200 ML Capacity, 1/4" I.D. Patient Tube, Pediatric Connectors, Sterile

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
<b>2012-320</b>	A8020L, A8020LF A8020-08L, A8020-08LF	N/A	N/A	6	<b>Baby Drain - ATS Compatible</b> ♦ Suction Control Stopcock ♦ In-Line Patient Tube Connector & Clamp

### Blood Recovery

♦ Large Capacity ATS Chamber, 1 Patient Tube, Sterile

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
<b>2050-000</b>	A9350L, A9350LF, A9350-08L, A9350-08LF	8888571364	N/A	6	<b>Blood Recovery - ATS Chamber</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock
<b>2050-000L</b>	A9350	N/A	N/A	6	<b>Blood Recovery - ATS Chamber</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock ♦ Latex Patient Tube
<b>2050-070</b>	N/A	N/A	N/A	6	<b>Blood Recovery - ATS Chamber</b> ♦ In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock ♦ Sterile Fluid Path Packaging



## Ocean Water Seal Chest Drains

### Blood Recovery

♦ Large Capacity ATS Chamber, 1 Patient Tube, Sterile

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
2050-300	A9350L, A9350LF, A9350-08L, A9350- 08LF	8888571364	N/A	6	<b>Blood Recovery - Standard System</b> ♦ No In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock

### Dual Blood Recovery

♦ Dual ATS Chamber, 2 Patient Tube, Sterile

Atrium Code #	Teleflex Pleur-evac	Tyco Aqua-Seal	Tyco Thora-Seal III	Case Qty.	Description
2052-000	N/A	N/A	N/A	6	<b>Dual Blood Recovery - ATS Chamber</b> ♦ 2 In-Line Patient Tube Connectors & Clamps ♦ Suction Control Stopcock
2052-300	N/A	N/A	N/A	6	<b>Dual Blood Recovery - ATS Chamber</b> ♦ No In-Line Patient Tube Connector & Clamp ♦ Suction Control Stopcock

## Oasis Dry Suction Control Chest Drains

### Single Collection Dry Suction Control

♦ 1 Patient Tube, Dry Suction Regulator, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
3600-100	A6000A, A6000L, A6000LF, A6000-08L, A6000-08LF	S1100C, S1100-08L, S1100-08LF	888571370	8888571489, 8888571562	6	<b>Single Collection - ATS Compatible</b> ♦ Pre-packaged water ♦ In-Line Patient Tube Connector & Clamp
3600-100L					6	<b>Single Collection - ATS Compatible</b> ♦ Pre-packaged water ♦ In-Line Patient Tube Connector & Clamp ♦ Latex Patient Tube
3600-150	A6050A, A6050L, A6050LF, A6050-08L, A6050-08LF	S1200, S1200-08L, S1200-08LF	888571370	8888571463	3	<b>Single Collection - ATS Ready</b> ♦ Pre-packaged water ♦ Preattached with In-Line ATS Blood Bag
3600-170	N/A	S1100FS, S1100-08FS	N/A	N/A	6	<b>Single Collection - ATS Compatible</b> ♦ Pre-packaged water ♦ In-Line Patient Tube Connector & Clamp ♦ Sterile Fluid Path Packaging

### Dual Collection Dry Suction Control

♦ 2 Patient Tubes, Major/Minor Collection Chambers, Dry Suction Regulator, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
3620-100	A6002L, A6002LF, A6002-08L, A6002-08LF	N/A	N/A	8888571513	6	<b>Dual Collection - ATS Compatible</b> ♦ Pre-packaged water ♦ 2 In-Line Patient Tube Connectors & Clamps

## Oasis Dry Suction Control Chest Drains

### Infant / Pediatric Dry Suction Control

◆ 200 ML Capacity, 1/4" I.D. Patient Tube, Pediatric Graduations and Connectors, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
3612-400	A6020L, A6020LF, A6020-08L, A6020- 08LF	N/A	N/A	N/A	6	<b>Baby Drain - Standard System</b> ◆ Pre-packaged water ◆ No In-Line Patient Tube Connector & Clamp
3612-100	A6020L, A6020LF, A6020-08L, A6020- 08LF	N/A	N/A	N/A	6	<b>Baby Drain - ATS Compatible</b> ◆ Pre-packaged water ◆ In-Line Patient Tube Connector & Clamp

### Blood Recovery Dry Suction Control

◆ 1 Patient Tube, ATS Chamber, Dry Suction Regulator, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
3650-100	A9250L, A9250LF, A9250-08L, A9250- 08LF	S1150C, S1150- 08L, S1150-08LF	N/A	N/A	6	<b>Blood Recovery - ATS Chamber</b> ◆ Pre-packaged water ◆ In-Line Patient Tube Connector & Clamp

### Dual Blood Recovery Dry Suction Control

◆ 2 Patient Tubes, Dual ATS Chamber, Dry Suction Regulator, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
3652-100	N/A	N/A	N/A	N/A	6	<b>Dual Blood Recovery - ATS Chamber</b> ◆ Pre-packaged water ◆ 2 In-Line Patient Tube Connectors & Clamps



## Express Dry Control Chest Drains

### Single Collection Dry Control

◆ 1 Patient Tube, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
4000-100N	N/A	S1100C, S1100-08L, S1100-08LF	N/A	N/A	6	<b>Single Collection - ATS Compatible</b> ◆ Pre-packaged Water ◆ In-Line Patient Tube Connector & Clamp
4000-170N	N/A	S1100FS, S1100-08FS	N/A	N/A	6	<b>Single Collection - ATS Compatible</b> ◆ Pre-packaged Water ◆ In-Line Patient Tube Connector & Clamp ◆ Sterile Fluid Path Packaging

### Blood Recovery Dry Suction Control

◆ 1 Patient Tube, ATS Chamber, Sterile

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Altitude	Tyco Sentinel Seal	Case Qty.	Description
4050-100N	N/A	S1150C, S1150-08L, S1150-08LF	N/A	N/A	6	<b>Blood Recovery - ATS Chamber</b> ◆ Pre-packaged Water ◆ In-Line Patient Tube Connector & Clamp
4050-200	N/A	S1150C, S1150-08L, S1150-08LF	N/A	N/A	6	<b>Blood Recovery - ATS Chamber</b> ◆ Pre-packaged Water ◆ Needleless Access Sampling Port ◆ In-Line Patient Tube Connector & Clamp
4050-170N	N/A	S1150C, S1150-08L, S1150-08LF	N/A	N/A	6	<b>Blood Recovery - ATS Chamber</b> ◆ Pre-packaged Water ◆ In-Line Patient Tube Connector & Clamp ◆ Sterile Fluid Path Packaging



**Pneumostat Chest Drain Valve**

♦ Sterile

Atrium Code #	B D Heimlich Valve	UreSil Tru Close	Arrow	Case Qty.	Description
16100	373460	TV11-10	AI - 01500E	10	<b>Chest Drain Valve</b> ♦ Stepped Connector w/Tube for Smaller Catheters ♦ One Way Valve ♦ 30ml Collection Volume ♦ Air Leak Well

**Express Mini 500 - Dry Suction , Dry Seal Mobile Chest Drain**

♦ Sterile

Atrium Code #	Case Qty.	Description
16400	6	<b>Mobile Chest Drain</b> ♦ No Water Required ♦ 500ml Capacity ♦ One Way Valve ♦ Connector Options, Straps
19930	10	Express Suction Bulb



## ATS Blood Bags and Accessories

### In-Line ATS Blood Bag

- ◆ For Use With Atrium ATS Compatible Chest Drains
- ◆ Sterile, Single Use
- ◆ 600 ML Capacity

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Aqua-Seal	Tyco Seal III	Thora- Seal	Tyco Sentinal Seal	Tyco Altitude	Case Qty.	Description
<b>2550</b>									<b>In-Line ATS Blood Bag</b>
Compatible With:	A1500	S100	8884713184	8884713184	8884713184	8884713184	8884713184	6	For use with the following ATS Compatible, In-Line Connector Models:
2002-000									Single Collection System
2002-100									Single Collection System
2002-050									Single Collection System
2020-000									Dual Collection System
2020-100									Dual Collection System
2050-000									Blood Recovery System
2050-100									Blood Recovery System
2052-000									Dual Blood Recovery System
3600-100									Dry Single Collection System
3600-150									Dry Single Collection System
3620-100									Dry Dual Collection System
3650-100									Dry Blood Recovery System
3652-100									Dry Blood Recovery System
4000-100N									Single Collection Dry Control
4000-170N									Single Collection Dry Control
4050-100N									Blood Recovery Dry Control System
4050-170N									Blood Recovery Dry Control System



## ATS Blood Bags and Accessories

### Self-Filling ATS Blood Bag

- ◆ For Use with Atrium Blood Recovery Chest Drains
- ◆ Sterile, Single Use
- ◆ 700 ML Capacity

Atrium Code #	Teleflex Pleur-evac	Teleflex Sahara	Tyco Aqua-Seal	Tyco Seal III	Thora-	Tyco Sentinal Seal	Tyco Altitude	Case Qty.	Description
<b>2450</b>									<b>Self Filling ATS Blood Bag</b>
Compatible With:	A1650	S250	8884571350	N/A		N/A	N/A	10	For use with the following ATS Chamber Models:
2050-000									Blood Recovery System
2050-100									Blood Recovery System
2050-300									Blood Recovery System
2052-000									Dual Blood Recovery System
2052-300									Dual Blood Recovery System
3650-100									Dry Blood Recovery System
3652-100									Dry Dual Blood Recovery System
4050-100N									Blood Recovery Dry Control System
4050-170N									Blood Recovery Dry Control System



## Chest Drain Accessories

Atrium Code #	Case Qty.	Description
19911	1	Standard Drain Caddy
19902	1	Large Drain Caddy
19907	10	Latex Free Tube Pack
19908	10	Latex Tube Pack
19909	10	Y - Extension Sets
19910	25	Red Cap Sterile
19913	25	Adult/Pediatric Connector 1/4" x 3/8"
19914	25	Standard Adult Connector 3/8" x 3/8"
19915	25	Pediatric Connector 1/4" x 1/4"
19916	25	Sims Connector 3/16" x 1/4"
19917	25	Y - Connector 3/8" x 3/8" x 3/8"
19918	25	Step - Down Connector 3/16" x 3/8"
19921	10	Needleless Access Sampling Port W/Tube
19922	100	Red Cap Bulk Non-Sterile
19923	100	Bulk Tubing - Latex
19924	100	Bulk Tubing - Non Latex
19925	100	Bulk Pack, Removable Clamp
19926	6	Dual Y Patient Line
19927	25	Connector Luer-Lock with Tube
19928	10	Y Suction Adapter
19932	20	Latex Free Tube Pack 1/4" x 4' (No Connector)



## PVC Taper-Tip Thoracic Catheters

♦ Sterile

♦ Thermosensitive PVC

♦ **Straight Catheter** - Tapered Connector Tip

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
8008	N/A	N/A	N/A	10	8	2
8010	N/A	N/A	N/A	10	10	2
8012	N/A	DSTC-12S	N/A	10	12	4
8016	N/A	DSTC-16S	500016	10	16	4
8020	N/A	DSTC-20S	500020	10	20	4
8024	N/A	DSTC-24S	500024	10	24	5
8028	8888570341	DSTC-28S	500028	10	28	6
8032	8888570358	DSTC-32S	500032	10	32	6
8036	8888570366	DSTC-36S	500036	10	36	6
8040	N/A	DSTC-40S	500040	10	40	6

♦ **Right Angle Catheter** - Tapered Connector Tip

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
8112	N/A	N/A	N/A	10	12	4
8116	8888570408	DRAC-16S	511016	10	16	4
8120	8888570416	DRAC-20S	N/A	10	20	4
8124	N/A	DRAC-24S	511124	10	24	4
8128	8888570432	DRAC-28S	511128	10	28	5
8132	8888570440	DRAC-32S	511132	10	32	5
8136	8888570457	DRAC-36S	511136	10	36	5
8140	8888570465	N/A	N/A	10	40	5



**ATP Trocar PVC Catheter - Includes Tapered Stylet**

<b>Atrium Code #</b>	<b>Kendall</b>	<b>Teleflex</b>	<b>Axiom</b>	<b>Case Qty.</b>	<b>Size (Fr)</b>	<b># Eyelets</b>
<b>8408</b>	8888560805	N/A	522208	10	8	2
<b>8410</b>	8888561019	DRTC-10S	522210	10	10	2
<b>8412</b>	8888561027	DRTC-12S	522212	10	12	2
<b>8416</b>	8888561035	DRTC-16S	522216	10	16	2
<b>8420</b>	8888561043	DRTC-20S	522220	10	20	2
<b>8424</b>	8888561050	DRTC-24S	522224	10	24	2
<b>8428</b>	8888561068	DRTC-28S	522228	10	28	2
<b>8432</b>	8888561076	DRTC-32S	855532	10	32	2

**HydraGlide™ Heparin Coated PVC Thoracic Catheters**

♦ **Straight Catheter - Tapered Connector Tip**

<b>Atrium Code #</b>	<b>Kendall</b>	<b>Teleflex</b>	<b>Axiom</b>	<b>Case Qty.</b>	<b>Size (Fr)</b>	<b># Eyelets</b>
<b>9008</b>	N/A	N/A	N/A	10	8	2
<b>9010</b>	N/A	N/A	N/A	10	10	2
<b>9012</b>	574202	N/A	N/A	10	12	4
<b>9016</b>	574210	N/A	N/A	10	16	4
<b>9020</b>	574228	N/A	N/A	10	20	4
<b>9024</b>	574236	N/A	N/A	10	24	5
<b>9028</b>	574244	N/A	N/A	10	28	6
<b>9032</b>	574251	N/A	N/A	10	32	6
<b>9036</b>	574269	N/A	N/A	10	36	6
<b>9040</b>	574277	N/A	N/A	10	40	6



## HydraGlide™ Heparin Coated PVC Thoracic Catheters

### ♦ Right Angle Catheter - Tapered Connector Tip

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
9112	N/A	N/A	N/A	10	12	4
9116	574301	N/A	N/A	10	16	4
9120	574302	N/A	N/A	10	20	4
9124	574327	N/A	N/A	10	24	4
9128	574335	N/A	N/A	10	28	5
9132	574343	N/A	N/A	10	32	5
9136	574350	N/A	N/A	10	36	5
9140				10	40	5

### ♦ ATP HydraGlide™ Trocar - Includes Tapered ATP Stylet

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
9408	N/A	N/A	N/A	10	8	2
9410	N/A	N/A	N/A	10	10	2
9412	N/A	N/A	N/A	10	12	2
9416	N/A	N/A	N/A	10	16	2
9420	N/A	N/A	N/A	10	20	2
9424	N/A	N/A	N/A	10	24	2
9428	N/A	N/A	N/A	10	28	2
9432	N/A	N/A	N/A	10	32	2

## Firm PVC Thoracic Catheters

- ◆ Open Tip for Fast, Secure Connection
- ◆ Thermosensitive PVC
- ◆ Sterile
- ◆ **Straight Catheter**

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
<b>15012</b>	8888570507	DSTC-12H	N/A	10	12	4
<b>15016</b>	8888570515	DSTC-16H	500016	10	16	4
<b>15020</b>	8888570523	DSTC-20H	500020	10	20	4
<b>15024</b>	8888570531	DSTC-24H	500024	10	24	5
<b>15028</b>	8888570549	DSTC-28H	500028	10	28	6
<b>15032</b>	8888570556	DSTC-32H	500032	10	32	6
<b>15036</b>	8888570564	DSTC-36H	500036	10	36	6
<b>15040</b>	8888570572	DSTC-40H	500040	10	40	6

### ◆ Right Angle Catheter Firm PVC Thoracic Catheter

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
<b>15212</b>	N/A	N/A	N/A	10	12	4
<b>15216</b>	8888571018	DRAC-16H	510016	10	16	4
<b>15220</b>	8888571026	DRAC-20H	N/A	10	20	4
<b>15224</b>	8888571034	DRAC-24H	511124	10	24	4
<b>15228</b>	8888571042	DRAC-28H	511128	10	28	5
<b>15232</b>	8888571059	DRAC-32H	511132	10	32	5
<b>15236</b>	8888571067	DRAC-36H	511136	10	36	5
<b>15240</b>	8888571075	DRAC-40H	511140	10	40	5



## Silicone Thoracic Catheters

- ♦ Sterile
- ♦ High Tear Strength Silicone
- ♦ Straight Catheter

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
14016	N/A	N/A	100016	10	16	4
14020	8888572537	N/A	100020	10	20	5
14024	8888572545	N/A	100024	10	24	6
14028	8888572552	N/A	100028	10	28	6
14032	8888572560	N/A	100032	10	32	6
14036	8888572578	N/A	100036	10	36	6

- ♦ Right Angle Catheter

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
14216	N/A	N/A	111016	10	16	4
14220	N/A	N/A	111120	10	20	5
14224	N/A	N/A	111124	10	24	6
14228	8888573048	N/A	111128	10	28	6
14232	8888573055	N/A	111132	10	32	6
14236	8888573063	N/A	111136	10	36	6



## Silicone Thoracic Catheters

- ◆ *Mediastinal Catheter*
- ◆ *Sterile*
- ◆ *High Tear Strength Silicone*

<i>Atrium Code #</i>	<i>Kendall</i>	<i>Teleflex</i>	<i>Axiom</i>	<i>Case Qty.</i>	<i>Size</i>	<i># Eyelets</i>
<b>14407</b>	N/A	N/A	111079	10	7mm	6
<b>14409</b>	N/A	N/A	111099	10	9mm	9
<b>14411</b>	8888573568	N/A		10	11mm	9
<b>14607</b>	N/A	N/A	111076	10	7mm	4
<b>14609</b>	N/A	N/A	111096	10	9mm	6
<b>14807</b>	N/A	N/A	111075	10	7mm	5
<b>14809</b>	8888573550	N/A	111095	10	9mm	5



**HydraGlide™ Heparin Coated Silicone Thoracic Catheters**

♦ *Straight HydraGlide™ Catheter*

♦ *Sterile*

♦ *High Tear Strength Silicone*

<b>Atrium Code #</b>	<i>Kendall</i>	<i>Teleflex</i>	<i>Axiom</i>	<b>Case Qty.</b>	<b>Size (Fr)</b>	<b># Eyelets</b>
<b>14116</b>	N/A	N/A	100016AT	10	16	4
<b>14120</b>	N/A	N/A	100020AT	10	20	5
<b>14124</b>	N/A	N/A	100024AT	10	24	6
<b>14128</b>	N/A	N/A	100028AT	10	28	6
<b>14132</b>	N/A	N/A	100032AT	10	32	6
<b>14136</b>	N/A	N/A	100036AT	10	36	6



**HydraGlide™ Heparin Coated Silicone Thoracic Catheters**

♦ **Mediastinal HydraGlide™ Catheter**

♦ Sterile

♦ High Tear Strength Silicone

♦ **Eyes on Top & Bottom**

<b>Atrium Code #</b>	<i>Kendall</i>	<i>Teleflex</i>	<i>Axiom</i>	<b>Case Qty.</b>	<b>Size</b>	<b># Eyelets</b>
<b>14507</b>	N/A	N/A	111079AT	10	7mm	6
<b>14509</b>	N/A	N/A	111099AT	10	9mm	9

♦ **Eyes on Top**

<b>Atrium Code #</b>	<i>Kendall</i>	<i>Teleflex</i>	<i>Axiom</i>	<b>Case Qty.</b>	<b>Size</b>	<b># Eyelets</b>
<b>14707</b>	N/A	N/A	111076AT	10	7mm	4
<b>14709</b>	N/A	N/A	111096AT	10	9mm	6

♦ **Eyes on Sides**

<b>Atrium Code #</b>	<i>Kendall</i>	<i>Teleflex</i>	<i>Axiom</i>	<b>Case Qty.</b>	<b>Size</b>	<b># Eyelets</b>
<b>14907</b>	N/A	N/A	111075AT	10	7mm	5
<b>14909</b>	573550	N/A	111095AT	10	9mm	5



## HydraGlide XL™ Heparin Coated Silicone Thoracic Catheters

- ◆ HydraGlide™ Catheter
- ◆ Sterile
- ◆ High Tear Strength Silicone
- ◆ Increased Number of Eyelets
- ◆ Increased Length by 6"

Atrium Code #	Ethicon Blake Wound Drain	Bard Channel Drain	Cardinal Hemaduct	Case Qty.	Size	# Eyelets
14720	2230 (19Fr)	072190 (19Fr)	CT-HU1900	10	20	15
14724	N/A	N/A	N/A	10	24	18

## PVC Taper-Tip Thoracic Catheters

- ◆ Thermosensitive PVC
- ◆ Straight Catheter - Tapered Connector Tip

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
8028-025	570549	DSTC-28S	500028	25	28	6
8032-025	570556	DSTC-32S	500032	25	32	6
8036-025	570564	DSTC-36S	500036	25	36	6
8040-025	570572	DSTC-40S	500040	25	40	6

- ◆ Right Angle Catheter - Tapered Connector Tip

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Size (Fr)	# Eyelets
8128-025	571042	DRAC-28S	511128	25	28	5
8132-025	571059	DRAC-32S	511132	25	32	5
8136-025	571067	DRAC-36S	511136	25	36	5



## Catheter Accessories

### ♦ PleuraGuide Disposable Chest Tube Kits

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Description
17100	565507	N/A	N/A	6	<b>Disposable Chest Tube Kit</b> ♦ All-in-one kit offers convenience and time savings ♦ Disposable - Eliminates cleaning and resterilization ♦ Collect kit along with preferred chest tube & drainage device
17150	565507	N/A	N/A	6	<b>Chest Tube Insertion Kit</b> ♦ No Lidocaine ♦ No Chloraprep ♦ CE Approved ♦ No # 10 scalpel
17200	565507	N/A	N/A	6	<b>PleuraGuide Kit</b> ♦ No Lidocaine

### ♦ PleuraGuide Chest Tube Introducer

Atrium Code #	Kendall	Teleflex	Axiom	Case Qty.	Description
19929	N/A	N/A	N/A	20	PleuraGuide Chest Tube Introducer